

98027

# Colour Coding Rubber Band x 5, Ø40 mm, Orange



Vikan rubber bands are the hygienic solution for secondary colour coding of cleaning equipment. The convenient bands fit shovels and vacuum cleaner nozzles.

# Technical Data

|  |                 |
|--|-----------------|
| <b>Item Number</b>   | 98027           |
| <b>Material</b>  | Silicone rubber |
| <b>Complies with (EC) 1935/2004 on food contact materials<sup>1</sup></b>              | No              |
| <b>Produced according to EU Regulation 2023/2006/EC of Good Manufacturing Practice</b> | Yes             |
| <b>FDA compliant raw material (CFR 21)</b>   | Yes             |
| <b>Use of phthalates and bisphenol A</b>   | No              |
| <b>Box Quantity</b>  | 1 Pcs.          |
| <b>Quantity per Pallet (80 x 120 x 200 cm)</b>   | 14000 Pcs       |
| <b>Quantity Per Layer (Pallet)</b>   | 1 Pcs.          |
| <b>Box Length</b>  | 100 mm          |
| <b>Box Width</b>   | 60 mm           |
| <b>Box Height</b>  | 40 mm           |
| <b>Product Diameter</b>  | 40 mm           |
| <b>Gross Weight</b>  | 0,02 kg         |
| <b>Net Weight</b>  | 0,02 kg         |
| <b>Cubik metre</b>   | 0 M3            |
| <b>Recommended sterilisation temperature (Autoclave)</b>                               | 121 °C          |
| <b>Max. cleaning temperature (Dishwasher)</b>  | 93 °C           |
| <b>Max usage temperature (food contact)</b>  | 100 °C          |
| <b>Max usage temperature (non food contact)</b>  | 200 °C          |
| <b>Min. usage temperature<sup>3</sup></b>  | -20 °C          |
| <b>Min. pH-value in usage concentration</b>  | 2 pH            |
| <b>Max. pH-value in Usage Concentration</b>  | 10,5 pH         |
| <b>Gtin-13 Number</b>  | 5705022002548   |
| <b>Customs Tariff No.</b>  | 40070000        |
| <b>Country of origin</b>   | China           |

New equipment should be cleaned, disinfected, sterilised and any labels removed, as appropriate to its intended use, e.g. high risk vs. low risk food production areas, general hospital areas vs. intensive care units, before use.