
Hygienic Wall Bracket, Grip Band Module, 83 mm, Blue





The Grip band module is designed for hanging cleaning tools without a hanging hole. Slide the Grip band module onto the supplied double-base/spacer from the left or right side. The Grip module can hold products with a diameter of 28-33 mm. The Grip band module is easy to disassemble for cleaning or replacement.

Technical Data

| Item Number | 10133 |
|---|--|
| Holds products between | 25 - 34 mm |
| Material | Polypropylene TPE Rubber Polyamide |
| Complies with (EC) 1935/2004 on food contact materials ¹ | No |
| Produced according to EU Regulation 2023/2006/EC of Good Manufacturing Practice | Yes |
| FDA compliant raw material (CFR 21) | Yes |
| Meets the REACH Regulation (EC) No. 1907/2006 | Yes |
| Use of phthalates and bisphenol A | No |
| Is Halal and Kosher compliant | Yes |
| Box Quantity | 5 Pcs. |
| Quantity Per Layer (Pallet) | 75 Pcs. |
| Box Length | 258 mm |
| Box Width | 230 mm |
| Box Height | 85 mm |
| Length | 83 mm |
| Width | 87 mm |
| Height | 48 mm |
| Net Weight | 0,08201 kg |
| Weight bag | 0,004 kg |
| Weight cardboard | 0,0232 kg |
| Tare total | 0,0272 kg |
| Gross Weight | 0,11 kg |
| Cubik metre | 0,000347 M3 |
| Recommended sterilisation temperature (Autoclave) | 121 °C |
| Max. cleaning temperature (Dishwasher) | 93 °C |
| Max usage temperature (non food contact) | 80 °C |
| Min. usage temperature ³ | -10 °C |
| Min. pH-value in usage concentration | 2 pH |
| Max. pH-value in Usage Concentration | 10,5 pH |
| Recycling Symbol "5", Polypropylene (PP) | No |
| Gtin-13 Number | 5705022030008 |
| GTIN-14 Number (Box quantity) | 15705028029997 |
| Customs Tariff No. | 39241000 |
| Country of origin | Denmark |

2/23/2021 2/3

New equipment should be cleaned, disinfected, sterilised and any labels removed, as appropriate to its intended use, e.g. high risk vs. low risk food production areas, general hospital areas vs. intensive care units, before use.

3. Do not store the product below 0° Celsius.

2/23/2021 3/3